# SOFARCLEAN®

Soluzione fisiologica sterile - Bomboletta pressurizzata da 150 ml/Sterile physiological solution - Pressurized can of 150 ml  
Solution physiologique stérile -Bouteille pressurisée de 150 ml/Sterile physiologische Kochsalzlösung - Druckluftspray mit 150 ml

## IT

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
## EN

### What you should know about SOFARCLEAN®

Read this product information leaflet carefully, before you use this product. It contains important information.

- SOFARCLEAN® can be purchased in pharmacy without prescription.
- The use of SOFARCLEAN® is very easy, but it must be used correctly to obtain the best results.
- Keep this leaflet. You may need to read it again.
- If you have any doubts or you need further information, contact the Medical Directorate of SOFAR, your doctor or pharmacist.

This leaflet explains the following information:

1. What is SOFARCLEAN®
2. Ingredients.
3. Operation.
4. Purpose.
5. Function.
6. When and how to use SOFARCLEAN®
7. Storage.
8. SOFARCLEAN® and the environment.
9. Other information.


#### 1. What is SOFARCLEAN®

SOFARCLEAN® is a medical device, consisting of a pressurised can containing sterile physiological solution.

#### 2. Ingredients

The can contains 150ml of sterile isotonic physiological solution, more specifically sterile distilled water containing 0.9% of sodium chloride, pressurised with compressed air at 9 bar. It does not contain any preservatives or propellant gases which are harmful for the ozone.

Content of the box: One 150ml can, 1 strong jet nozzle (blue colour), 1 vaporised jet nozzle (white), product information leaflet.

#### 3. Operation

SOFARCLEAN® is a pressurised can which uses a new system to force-supply sterile physiological solution. The innovative quality of this system is that the propellant is only compressed air, which is totally separated from the water solution contained in the can. The forced supply of the solution is obtained by simply pressing the nozzle of the can (the application is shown in the drawing). Pressure on the nozzle produces a continuous flow, in whatever position (at 360°) the can is held. The flow is interrupted by taking the pressure off the nozzle.

#### 4. Purpose

SOFARCLEAN® has been designed to CLEANSE, IRRIGATE and HYDRATE acute and chronic wounds, dermatological lesions of various kinds, cavities and mucosa of any body part, healthy or injured cutaneous surfaces.

#### 5. Function

SOFARCLEAN® comes with two different nozzles, each with its own function, described in point A and B:

A) The first blue-coloured nozzle produces a strong compact jet, with an elevated mechanical cleansing effect, particularly useful for removing foreign bodies, dirt, necrotic tissue, serum exudate, etc. from the wound. The strong jet of SOFARCLEAN® is also useful for cleansing acute wounds, essential in first aid, also as an aid in the mechanical removal of crusty necrotic tissue. The power of the jet is adjusted by bringing the "beak" closer or further away from the wound, to reach the most appropriate distance for the purpose. The same jet can also be used as an irrigator, for washing those parts of the body which have become inaccessible to daily hygiene, due to transitory or permanent immobility of the patient due to accident or paralysis.

B) The second white-coloured nozzle produces a vaporised jet, suitable for gently hydrating, irrigating or cleansing: mucosa, noses, throat, mouth, ear, anogenital area, surface injuries, dermatological lesions of various kind. SOFARCLEAN® is also safe and harmless and can be used in individuals of any age, including small children, without risk. The solution is supplied at ambient temperature, so it does not cool down the surface to which it is applied (phenomenon connected with fast evaporation propellants). The use of SOFARCLEAN® excludes any type of allergic reaction, as it is free of preservatives, impurities, foreign substances, propellant gases.

#### 6. When and how to use SOFARCLEAN®

In the local treatment of wounds, the key measures to be taken are two: first cleanse adequately, and if necessary, disinfect appropriately, two procedures which should be considered separately. It is wrong to see the two, cleansing and disinfection, as a combined action, using the many products available on the market which claim to have this double function. SOFARCLEAN® adequately responds to the essential primary need of cleansing acute and chronic wounds and dermatological lesions of various kind. Its application is virtually harmless and without any undesired side-effects. As regards the second measure, disinfection of the wound, it must be considered that it can very often cause chemical traumas in the wound itself, seriously compromising the repair process. By definition, a wound is a lesion which must heal as quickly as possible. Further traumas, caused involuntarily in the attempt to medicate it, will inevitably delay the healing time. Traumas which delay or sometimes prevent the healing of an wound, can be divided into two large categories: mechanical and chemical traumas. Mechanical traumas, include cleansing which is too fast, careless rubbing on the bed of the wound or careless removal of the adherent bandage, etc. Chemical traumas inflicted on wounds are much greater than the mechanical ones, because caused by the indiscriminate use of disinfectants, above all in liquid form, used with the double effect of cleanser and disinfectant. All these disinfectants in liquid form, none excluded, are physiologically incompatible with the wound, they are sometimes caustic, and are in every case indistinctly cytotoxic. The components of the various products available on the market are toxic for the tissue defence systems, harmful for the wound, and inhibit its healing. An American research professor has believed these elementary concepts for years, corroborated by experimental tests: G. Rodeheaver; *Controversies in the topical treatment of wounds*; vol. 1 (1) April 1989 pg. 19-27. In his own words: "Disinfectants of any kind damage the tissues and interfere with the tissue function, aggravating the wound and delaying its healing." In the topical treatment of wounds, Rodeheaver recommends only the use of sterile physiological solution, decisively condemning, simply because they are protoplasmic killers, all disinfectants based on iodine, iodopovidone, oxygenated water, chlorhexidine, boric acid, various alcohols, exachlorophene, formaldehyde, hypochlorite, acetic acid, Silver nitrate, merthiolate, gentian violet, permanganate, aluminium salts. This is also recommended by the American guidelines on the treatment of lesions from bedsore. Only when the wound does not tend to heal, or there are clear signs of infection should a medicine for topical antibacterial use be applied.

**SOFARCLEAN® with blue nozzle, strong jet, is recommended for the systemic cleaning of acute and chronic wounds, before applying the medication. A careful cleaning, by irrigation with SOFARCLEAN® allows to appropriately medicate the lesions, without resorting to disinfectants. The application of SOFARCLEAN® does not produce any cold feeling on the wound, because it is at ambient temperature and free of rapid evaporation propellants.**

**IMPORTANT:** before applying the jet to the wound, always spray a little solution elsewhere in order to cleanse the nozzle.

**SOFARCLEAN® with white nozzle, vaporised jet, is suitable for gently hydrating, irrigating or cleansing: nostrils, throat, mouth, ear, anogenital areas, mild and surface injuries, dermatological lesions of various kinds. In the event of illness from cold, SOFARCLEAN® is particularly effective for freeing the nasal cavities, in the adult, child and infant, from the presence of obstructive mucus, which often disturbs breathing and sleep. Nasal irrigation, recommended at least four times a day, contributes to freeing the nasal cavities, and soothing inflammation of the mucosa, to restore physiological respiration and reduce the local bacterial load to a minimum.**

**IMPORTANT:** before applying the vaporised jet, spray a little of the solution elsewhere to clean the nozzle.

#### 7. Storage

The can comes with a cap, which must be put back after use, leaving the nozzle in place. The second nozzle should be stored beside the can, in a closed box. The product does not require any specific storage, simply keep it at ambient temperature and out of children's reach. Before a new application, always remember to spray a little of the solution elsewhere, to clean the selected nozzle.

#### 8. SOFARCLEAN® and the environment

The components of SOFARCLEAN®, sterile physiological solution and air, are 100% eco-friendly and therefore respect the ecological balance.

#### 9. Other information

SOFARCLEAN® is proposed as an essential tool to introduce a highly innovative qualified procedure in the treatment of acute and chronic wounds able to determine the abandonment of the traditional and indiscriminate use of so-called disinfectants in general, which aside from having little effect, being cytotoxic, inevitably damage the wound repair processes. To prefer the sterile physiological solution of SOFARCLEAN® to disinfectants in medicating injuries, means assuring a physiological repair process, the absence of undesired reactions, lower costs and a faster acceleration of healing times.

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SOFAR S.p.A.