SOFARGEL®

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Silver Sulfadiazine GEL-based dressing

Product description: SofarGEL is a water-based dressing for topical use containing Silver Sulfadiazine. The high water content (over 95%) GEL keeps the wound hydrated, promoting the tissue healing process. It creates a mechanical barrier against external agents and protects the damaged skin due to the silver-based compound. It also promotes debridement in the event of tissue necrosis. The presence of Silver Sulfadiazine, a metal-organic compound obtained by the reaction of silver nitrate with sulfadiazine, prevents microbial contamination of the dressing and keeps the wound environment free from exogenous bacteria. The special formulation of SofarGEL creates an environment that promotes the natural healing process and provides an effective barrier to microbial penetration. The product remains easily in place and is not affected by the pulling of the wound edges. It can be easily removed by cleaning the wound with saline solution.

Composition: Purified Water, Glycerol, Carbopol 974P, Sodium Hydroxide, Micronized Silver Sulfadiazine

Indications: SofarGEL is indicated for the local treatment of skin wounds such as: abrasions, excoriations, cuts, wounds, skin irritation and dermatological conditions in general. It can be used as a moisturising agent to keep the wound moist and to prevent infections or as a cleansing agent to remove dead tissue and fibrin deposits.

Dosage and directions for use: After cleansing the wound properly, apply an even layer of SofarGEL, about 2 or 3 mm thick, on the affected surface. The GEL can be applied directly onto the wounds or previously spread on a sterile lint to achieve the moisturising effect or on a polyurethane or semioclusive dressing, when used as a cleansing agent. SofarGEL is also an efficient substitute for most commonly used *hydrogels*. The GEL must be applied continuously once or twice a day as long as the infection risk persists and until the wound is completely healed, both in the case of natural healing and after surgical interventions. Carefully clean the wound with water or saline solution everytime before reapplying the product.

Precautions: Patients with known hypersensitivity to its ingredients, must not use SofarGEL. It is not recommended to use the

product for more than 4 weeks. Consult a physician if the wound does not heal completely or if undesired symptoms occur. In the event of accidental ingestion, seek medical advice.

Warnings: The gel may tend to turn brown when exposed to light, due to the silver ions in the formulation. This occurrence, however, does not affect its efficacy and safety. Do not apply the product together with liquid or solid disinfectants. Close the tube after every use. The product should not be used in pregnant women or infants unless recommended by a physician. Sofargel is to be used for treating minor wounds only. The product is not intended for ophthalmic use.

KEEP THE PRODUCT OUT OF THE REACH AND SIGHT OF CHILDREN.

Adverse reactions: Temporary irritation or burn sensation may occur, although without consequences, and rare local allergic reactions. If any undesirable effect occurs, consult a physician or pharmacist.

Expiry date and storage: Store at temperature below 40 °C. Do not use the product after the expiry date indicated on the package. This date refers to the product properly stored in original packaging.

Packaging: Aluminium tube of 3g – 10g – 30g – 50g – 120g – 180g

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